

We're not just blowing hot air...



SCIP Approved “Air Free” Conductive Underbody Patient Warming Technology

- ▶ Reusable and energy efficient. Eliminates costly disposables, storage and waste.
- ▶ Warms every patient effectively and affordably – promoting normothermia.
- ▶ Eliminates concerns and inconveniences of blowing air/circulating water systems.



The Affordable Answer to Active Patient Warming...




As a specialist in critical care temperature monitoring, NOVAMED USA focuses on new technologies to minimize risk of hypothermia. Utilizing carbon polymer technology, KOALA^{ΔT} is the simple solution to warm every OR patient before, during, and after surgery.



- ▶ Actively warms to prevent perioperative hypothermia.
- ▶ Ensures complete patient access, ideal for open heart.
- ▶ X-Ray and radio translucent.
- ▶ Temperature settings available in:
 - °C: 37°, 38°, 39°, 40°
 - °F: 99°, 100°, 102°, 104°
- ▶ Operates in total silence. Noise-free.

Benefits of Koala^{ΔT} Warming System's Multi-layer Mattress



Pressure-Relief "Memory Foam"

- ▶ Conforms to improve patient circulation
- ▶ Provides even distribution of weight and warmth
- ▶ Prevents decubitus ulcers

Durable Latex-Free Cover

- ▶ Welded seams ensure long life
- ▶ Maintenance-Free
- ▶ Allows easy disinfection between patients

Comfort Layer

- ▶ Soothing layer of cotton provides enhanced comfort and pressure-relief

Patented Carbon Polymer Fabric

- ▶ Provides active warmth to effectively maintain normothermia
- ▶ Conductive thermal fabric provides continuous, evenly distributed warmth

Available Sizes

A blue circular logo with a diagonal slash through it, and the word "LATEX" in bold blue capital letters across the center.



Available to suit a full range of surgical procedures.

½ Length



Full Length



Available Sizes



Product Description	Part Number	Dimensions
KOALA Control Unit, 110V	29 - KOR - CON	
Standard Width Mattress Selection		
KOALA O.R. Mattress, Full Length, Narrow	29 - KOR - FLN	75" x 21"
KOALA O.R. Mattress, ¾ Length, Narrow	29 - KOR - MLN	47" x 21"
KOALA O.R. Mattress, ½ Length, Narrow	29 - KOR - HLN	42" x 21"
KOALA O.R. Mattress, Neonatal	29 - KOR - NEO	22" x 20"
Extended Width Mattress Selection		
KOALA O.R. Mattress, Full Length, Wide	29 - KOR - FL	75" x 23"
KOALA O.R. Mattress, ¾ Length, Wide	29 - KOR - ML	47" x 23"
KOALA O.R. Mattress, ½ Length, Wide	29 - KOR - HL	42" x 23"

NOVAMEDUSA

800-425-3535 • www.novamed-usa.com