

# NEOTEMP

**Neonatal Skin Temperature Sensors  
for the NICU, Newborn Nursery and L & D**



**For Precise Temperature Measurement  
with Baby Sensitive Accuracy**

Select from a complete range of disposable and reusable sensors designed to fit most models of GE/Ohmeda Medical and Dräger Medical/Air Shields Incubators & Infant Radiant Warmers.

**NOVAMEDUSA**

800-425-3535  
[www.NovaMed-USA.com](http://www.NovaMed-USA.com)

# NEOTEMP

**Precision Neonatal Temperature Monitoring  
for GE/Ohmeda Medical and Dräger Medical/Air Shields  
Incubators & Infant Radiant Warmers.**



- Temp Accuracy** Baby-sensitive thermistor sensor guarantees temperature reading sensitivity with a stable, highly accurate measurement.
- Compatibility** Choose from a comprehensive range of sensors designed for GE / Ohmeda Medical and Dräger Medical / Air Shields incubators and radiant warmers.
- Ease of Use** Molded plug assembly provides direct connection to the incubator/ warmer eliminating the need for an interface cable.
- Sensor Safety** Flexible lead wire facilitates sensor placement directly on the infant further ensuring precise temperature reading.
- Gentle** Ensures patient safety with a one inch round baby sensitive hypoallergenic adhesive cover. A reflective Mylar adhesive cover provides insulation from heat and external sources.

NEONATAL Replacement Adhesive Mylar Disc Covers		
10 - 1900 - 300	Mylar Adhesive Replacement Disc Covers -	300 / Box

# Precise Temperature Monitoring for Transport Incubators

**AIRBORNE**  
Life Support Systems

DM/AIRSHIELDS®  
500 TRANSPORT  
INCUBATOR



**400 Series  
Skin Temperature  
Sensors**

## Accuracy

NOVAMED USA's 400 Series thermistor sensor guarantees a continuous temperature reading with accuracy  $\pm 0.2^{\circ}\text{C}$  for the full range ( $25^{\circ}\text{C} - 50^{\circ}\text{C}$ ).

## Quick Response

Unique thermistor design ensures rapid response to reflect any temperature change on transport patients.

## Stress-Free

NOVATEMP sensors are engineered with a longer lead wire and a custom fitted molded connector to ensure response time reliability.

### NOVATEMP SKIN TEMPERATURE SENSORS / 400 SERIES

Order Ref #	Product	Packaging
10-1600-030	Skin Temperature Sensors/ 400 Series	30 / Box
10-1600-090	Skin Temperature Sensors/ 400 Series	90 / Case

## NOVATEMP Temperature Adapter Cable for AirBorne and Air Shields Transport Equipment



Reusable Temperature Adapter Cable with  $\frac{1}{4}$ " Right Angle Phono Plug Connector for AirBorne and Air Shields Transport Equipment and Molex Connector for 400 Series Temperature Sensors.

**Available in both 1' and 12' Lengths**

Order Ref # 11 - CA - 400 - 01 or #11- CA - 400 - 12

# NEOTEMP

## Disposable Neonatal Skin Temperature Sensors

(Packaged 20/box, 60/case)

		Part #	Description	Substitute For
GE/Ohmeda Medical/Ohio		10 - 1930 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for GE Care Plus, I.C., CG Incubators & Infant Warming Systems (IWS) 3000, 3300, 3500, 5000 & Panda Warmers	Replaces OEM # LA003, LA005, # 6600 - 0196-700 & # 6600 - 0208-700
		10 - 1930 - 060 Packaged 60/case		
GE/Ohmeda Medical/Ohio		10 - 1950 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for GE GIRAFFE OmniBed, Giraffe Incubators & Warmers, New Panda Warmers	Replaces OEM # 6600 - 0873-700 & # 6600 - 0874-700
		10 - 1950 - 060 Packaged 60/case		
Dräger Medical/Air-Shields		10 - 1910 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for ISOLETTE Incubators (Models C-100, C-450, C-550) & PM-78 Radiant Warmers	Replaces OEM #68-209-20/MU06933 & #68-209-30/MU06935
		10 - 1910 - 060 Packaged 60/case		
		10 - 1960 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for ISOLETTE C2000 & 8000 Incubators, VERSALET 7700 Care Center	Replaces OEM # 83-104-15/MU12551
		10 - 1960 - 060 Packaged 60/case		
	10 - 1965 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for CALEO Incubators	Replaces OEM # MX11000	
	10 - 1965 - 060 Packaged 60/case			
	10 - 1970 - 020 Packaged 20/box	NEOTEMP Neonatal Skin Temperature Sensor for RESUSCITAIRE Radiant Warmers	Replaces OEM # 81-300-006, # 81-300-008/MU11228 & # 81-300-009/MU11229	
	10 - 1970 - 060 Packaged 60/case			
Atom Medical		10 - 2000 - 020 Packaged 20/box	NOVATEMP Skin Temperature Sensor For Atom Medical Incubators & Infant Warmers	Replaces OEM #60882, #60883, & #60084 w / connector
	10 - 2000 - 060 Packaged 60/case			